

Alpha Thalassemia

Alpha (α) thalassemia is the most common inherited disorder of hemoglobin (Hb) worldwide and is caused by *HBA1* and *HBA2* gene variants. Decreased or absent synthesis of the hemoglobin (Hb) α chain results in clinical presentations ranging from asymptomatic silent carriers to severe anemia and fetal lethality. The two clinically significant forms of α thalassemia are Hb Bart hydrops fetalis syndrome and hemoglobin H (HbH) disease. Alpha thalassemia is found more often in certain ethnicities, including African, African American, Mediterranean, Middle Eastern, and Southeast Asian.

Disease Overview

Prevalence and/or Incidence

- Most common inherited disorder of Hb worldwide
- Carrier frequencies in high-risk populations:
 - African, African American: 1/3
 - Middle Eastern, Southeast Asian: 1/20
 - Mediterranean: 1/30-50
- Hb Bart hydrops fetalis syndrome and HbH disease are more frequent in Southeast Asian, Asian Indian, and Mediterranean populations than in African populations.

Symptoms

Phenotype	Associated Symptoms
α Thalassemia silent carrier	Typically asymptomatic, though borderline anemia or mild microcytosis may be present Often misdiagnosed as iron deficiency Normal Hb electrophoresis
α Thalassemia trait	Mild microcytic anemia may be present Often misdiagnosed as iron deficiency Normal Hb electrophoresis
HbH disease	Moderate to severe form of α thalassemia Moderate microcytic hypochromic anemia Hemolysis with Heinz bodies Splenomegaly Rare extramedullary hematopoiesis Propensity for acute hemolysis after oxidative stress, drug therapy, or infection

Tests to Consider

Alpha Thalassemia (*HBA1* and *HBA2*) Deletion/Duplication with reflex to Hb Constant Spring 3003651

Method: Multiplex Ligation-dependent Probe Amplification with Sanger sequencing confirmation of HbCS

- Preferred first-tier genetic test for confirmation of suspected α thalassemia or α thalassemia trait
- Detects common, rare, and novel deletions or duplications in the α -globin gene cluster and its HS-40 regulatory region
- Multiplex ligation-dependent probe amplification (MLPA) used to detect Hb Constant Spring (HbCS) (*HBA2* c.427T>C; p.Ter143Gln); targeted Sanger sequencing is performed to assess HbCS copy number in absence of a concurrent *HBA2* deletion.

Alpha Globin (*HBA1* and *HBA2*) Deletion/Duplication (Extended TAT as of 11/20/20-no referral available) 2011622

Method: Multiplex Ligation-dependent Probe Amplification

- First-tier genetic test for confirmation of suspected α thalassemia or α thalassemia trait
- Detects common, rare, and novel deletions or duplications in the α -globin gene cluster (*HBZ*, *HBM*, *HBA1*, *HBA2*, *HBQ1*) and its HS-40 regulatory region

Alpha Globin (*HBA1* and *HBA2*) Sequencing and Deletion/Duplication (Temporary Referral as of 05/10/21) 2011708

Method: Polymerase Chain Reaction/Sequencing./Multiplex Ligation-dependent Probe Amplification.

- Comprehensive genetic test for detection of α thalassemia or α thalassemia trait
- Detect deletional and nondeletional variants in *HBA1* and *HBA2*

Alpha Thalassemia (*HBA1* and *HBA2*) Deletion/Duplication with reflex to Hb Constant Spring, Fetal 3003656

Method: Multiplex Ligation-dependent Probe Amplification with Sanger sequencing confirmation of HbCS

- Diagnostic testing for α thalassemia in fetus with suggestive clinical findings or at risk for α thalassemia due to familial *HBA1*/*HBA2* deletions or HbCS variant

Phenotype	Associated Symptoms
Hb Bart hydrops fetalis syndrome	<p>Most severe form of α thalassemia</p> <p>Risk for fetus</p> <ul style="list-style-type: none"> Lethal in fetal or early neonatal period Generalized edema, ascites, pleural and pericardial effusions Severe hypochromic anemia Usually detected on ultrasound at 22-28 weeks gestation <p>Maternal complications during pregnancy</p> <ul style="list-style-type: none"> Preeclampsia Polyhydramnios or oligohydramnios Antepartum hemorrhage Premature delivery

- Detects common, rare, and novel deletions or duplications in the α -globin gene cluster and its HS-40 regulatory region
- MLPA is used to detect HbCS (*HBA2* c.427T>C; p.Ter143Gln); targeted Sanger sequencing is performed to assess HbCS copy number in absence of a concurrent *HBA2* deletion.

See [Related Tests](#)

Pathophysiology

Typically, individuals have four functioning α -globin genes ($\alpha\alpha/\alpha\alpha$). Two genes, *HBA1* and *HBA2*, are present on each copy of chromosome 16 and α -globin chains function as subunits of fetal Hb (HbF: $\alpha_2\gamma_2$) and adult Hb (HbA: $\alpha_2\beta_2$). The number of α -globin genes deleted or inactivated correlates with different α thalassemia phenotypes. Genotype/phenotype correlations in α thalassemia are complex and may be influenced by coinheritance of other Hb variants or α -globin gene duplications.

Phenotype	Genotype(s)
α Thalassemia silent carrier	$-\alpha/\alpha\alpha$
α Thalassemia trait	$-\alpha/-\alpha$ $--/\alpha\alpha$
HbH disease	$-\alpha/-$
Hb Bart hydrops fetalis syndrome	$--/--$

Genetics

Genes

HBA1 and *HBA2*

Inheritance

Autosomal recessive

Variants

- HBA1* and *HBA2* large gene deletions account for approximately 90% of pathogenic α -thalassemia variants.
 - $-\alpha 3.7$ and $-\alpha 4.2$ deletions result in the deletion of a single gene.
 - $-(\alpha)20.5$, $--SEA$, $--MED-I$, $--FIL$, and $--THAI$ deletions result in the deletion of both *HBA1* and *HBA2* genes from the same chromosome
- Sequence variants and regulatory region variants occur mainly in *HBA2* and account for up to 15% of causative variants.
 - Nondeletional variants include:
 - Sequence variants that inactivate the gene

- Small insertions/deletions
- Variants that result in unstable α -globin protein (eg, Hb Constant Spring or HbCS)
- Nondeletional α -globin variants may be pathogenic or benign.
 - Both may result in an abnormal protein detectable by Hb evaluation.
 - Pathogenic nondeletional variants often have a more severe effect than single gene deletions.
- Alpha-globin gene duplication results in three or more active α -globin genes on a single chromosome.
 - Typically benign
 - May alter expected clinical phenotypes and hematological features when coinherited with beta (β) thalassemia

Test Interpretation

Alpha Globin (*HBA1* and *HBA2*) Sequencing and Deletion/Duplication

Sensitivity/Specificity

- Analytical sensitivity/specificity: 99% for both duplication/deletion analysis and sequencing
- Clinical sensitivity: most pathogenic *HBA1* and/or *HBA2* gene variants are large deletions not detectable by sequencing
 - Deletion: at least 90%, varies by ethnicity¹
 - Sequencing: up to 15%, varies by ethnicity¹

Results

Result	Variant(s) Detected	Interpretation
Positive	Genotype (-a/aa)	Silent carrier
	Genotype (-a/-a) or (--/aa)	α Thalassemia trait
	1 pathogenic sequence variant	Silent carrier or α thalassemia trait
	Genotype (-a/aa)	HbH disease
	Genotype (--/-)	Hb Bart hydrops fetalis syndrome
	Genotype (aa/aaa)	Extra functional alpha globin gene copy is present
Negative	No pathogenic variants detected	Greatly decreased probability that the individual is affected with, or a carrier of, α thalassemia
Inconclusive	Variant(s) of uncertain clinical significance identified	Unknown clinical significance

Limitations

- Diagnostic errors can occur due to rare sequence variations.
- Sequence analysis will not detect all regulatory region variants or variants in α -globin cluster genes other than *HBA1* and *HBA2*.
- Sequencing of both *HBA1* and *HBA2* may not be possible in individuals harboring large α -globin deletions on both alleles.
- This assay is unable to sequence *HBA2-HBA1* fusion genes; thus, *HBA1* or *HBA2* sequence variants occurring in cis with a 3.7 kb deletion or other *HBA2-HBA1* hybrid gene will not be detected.
- It may not be possible to determine the phase of identified sequence variants.

- Specific breakpoints of large deletions/duplications will not be determined; therefore, it may not be possible to distinguish deletion variants of similar size.
- Individuals carrying both a deletion and duplication within the α -globin gene cluster may appear to have a normal number of α -globin gene copies.
- Rare syndromic or acquired forms of α thalassemia associated with *ATRX* variants will not be detected.

References

1. Origa R, Moi P. [Alpha-thalassemia](#). In: Adam MP, Ardinger HH, Pagon RA, et al, editors. GeneReviews, University of Washington; 1993-2021. [Last update: Oct 2020; Accessed: Jan 2021]

Related Tests

[Hemoglobin Evaluation with Reflex to Electrophoresis and/or RBC Solubility 0050610](#)

Method: High Performance Liquid Chromatography/Electrophoresis/RBC Solubility

[Hemoglobin Evaluation Reflexive Cascade 2005792](#)

Method: High Performance Liquid Chromatography/Electrophoresis/RBC Solubility/Polymerase Chain Reaction/Fluorescence Resonance Energy Transfer/Sequencing

[Familial Mutation, Targeted Sequencing 2001961](#)

Method: Polymerase Chain Reaction/Sequencing

[Familial Mutation, Targeted Sequencing, Fetal 2001980](#)

Method: Polymerase Chain Reaction/Sequencing

ARUP Laboratories is a nonprofit enterprise of the University of Utah and its Department of Pathology, 500 Chipeta Way, Salt Lake City, UT 84108
(800) 522-2787 | (801) 583-2787 | aruplab.com | arupconsult.com
Content Review January 2021 | Last Update May 2021