



HEMOSTASIS/THROMBOSIS SPECIMEN COLLECTION AND HANDLING

To produce valid results for hemostasis/thrombosis tests and factor assays, specimen integrity is crucial and must be maintained. All specimens sent for testing must be collected and shipped in the following manner:

Note: Anticoagulant medications often interfere with coagulation studies and this should be considered when ordering and interpreting coagulation tests.

1. Obtain venous blood by clean venipuncture at a site away from an intravenous line. Avoid slow-flowing draws and/or traumatic venipunctures, as either of these may result in an activated or clotted specimen. Do not use needles smaller than 23 gauge. Do not leave the tourniquet on for an extended length of time before drawing the specimen.
2. Draw a discard tube (nonadditive or light blue-top tube) as a recommended procedure before drawing coagulation specimens using a winged blood collection device in light blue-top vacuum tubes (3.2% buffered sodium citrate). The discard tube does not need to be completely filled; its purpose is to ensure that the collection device tubing is filled with blood. Discard or pilot tubes are not required for collection with other blood collection devices. Note: Reference intervals have been established using 3.2% buffered sodium citrate.
3. Fill light blue-top tubes as far as vacuum will allow and mix by gentle inversion five to six times. Exact ratio of nine parts blood to one part anticoagulant must be maintained. Inadequate filling of the specimen tube will alter this ratio and may lead to inaccurate results.

Note: Patients who have hematocrit values above 55% should have the anticoagulant adjusted to maintain the nine to one ratio by removing a portion of the volume of the citrate solution. To calculate the amount of citrate required in the collection tube or syringe, use the following formula:

$$C = (1.85 \times 10^{-3}) (100 - \text{HCT}) (V_{\text{Blood}})$$

Where:

- C (mL) is the volume of 3.2% citrate remaining in the tube.
- HCT (%) is the hematocrit of the patient.
- V is the volume of blood added (if a 2.7 mL tube is used, then the volume is 2.43 mL to achieve a nine to one blood to anticoagulant ratio).
- 1.85×10^{-3} is a constant.

Example calculation: Patient hematocrit is 60% and blood will be drawn into a 2.7 mL citrate tube.

Calculate the adjusted citrate volume:

$$C \text{ (mL)} = (1.85 \times 10^{-3}) (100 - 60) (2.43 \text{ mL}) = 0.18 \text{ mL}$$

Then calculate the amount of citrate to remove:

$$0.27 \text{ mL} - 0.18 \text{ mL} = 0.09 \text{ mL of citrate to remove.}$$

4. Centrifuge the specimen at 1700 x g for 15 minutes or at a speed and time required to consistently produce platelet-poor plasma (i.e., platelet count less than 10,000/ μ L). Hemolyzed specimens will be rejected.
5. Immediately remove only the top two-thirds of the platelet-poor plasma from the specimen using a plastic transfer pipette. (Use of glass transfer pipettes may result in activation and/or clotting of the plasma.) Place the plasma in a properly labeled ARUP standard transport tube and clearly mark the vial contents as plasma. Glass vials will be rejected.
6. Immediately freeze the plasma in a non-frost-free freezer. A frost-free freezer should not be used due to the warming cycles. Specimens may be stored at -20°C or -70°C. Refer to individual test information for frozen stabilities. Specimens must remain frozen during storage and shipment. Specimens should be stored capped. A separate ARUP standard transport tube must be submitted for each assay requested.

- Room temperature stabilities should be observed if a specimen cannot be frozen immediately.
- Uncentrifuged, centrifuged with plasma on top of cells, or centrifuged with plasma separated from cells should be kept at 18–24°C for:

Specimen Type	Room Temperature Stability
PT specimens	no longer than 24 hrs from the time of collection
aPTT specimens	no longer than 4 hrs from the time of collection

*not drawn to monitor heparin therapy

- Specimens for other coagulation tests have variable stabilities and should be stored in the same manner as aPTT specimens.
- Coagulation specimens should not be stored at refrigerated temperatures (2–8°C) or on ice.