

Acute Myeloid Leukemia Mutation Panel by Next Generation Sequencing

Last Literature Review: August 2023 Last Update: March 2026

Acute myeloid leukemia (AML) is a myeloid malignancy (i.e., a clonal disorder of hematopoietic stem and progenitor cells). Recent studies have identified recurrently mutated genes with diagnostic and/or prognostic significance in AML. The presence of certain variants may inform clinical management.

This multigene panel by massively parallel sequencing (also referred to as next generation sequencing [NGS]) detects molecular changes, including single nucleotide variants (SNVs) and small insertions and deletions (indels) relevant to AML. This panel is more cost-effective than using multiple single gene tests and can be used to complement the morphologic and cytogenetic workup of AML.

Disease Overview

Diagnostic, Prognostic, and Treatment Issues

Identification of one or more clonal genetic abnormalities, variants, or patterns of variants in patients with AML may:

- Aid in establishing the diagnosis and subclassification
- Aid in determining prognosis
- Inform clinical management

For more information about the recommended testing strategy for AML, refer to the ARUP Consult [Acute Myeloid Leukemia – AML](#) topic.

Genetics

Genes Tested

This panel tests *ANKRD26*, *ASXL1*, *CEBPA*, *DDX41*, *DNMT3A*, *ETV6*, *FLT3*, *GATA2*, *IDH1*, *IDH2*, *KIT*, *KRAS*, *NPM1*, *NRAS*, *RUNX1*, *TP53*, and *WT1*.

For some genes, one or more exons of the preferred transcript are not covered by sequencing for the indicated gene. Refer to the [Genes Tested](#) table below for a full list of targeted regions and exclusions.

Test Interpretation

Results

Reported variants are classified into two categories:

- Tier 1: Mutations with known clinical significance in hematologic malignancies
- Tier 2: Variants of unknown clinical significance in hematologic malignancies

Featured ARUP Testing

[Acute Myeloid Leukemia Mutation Panel by Next Generation Sequencing 3002714](#)

Method: Massively Parallel Sequencing

For more information on ARUP's myeloid malignancies panel, which tests the genes in this panel and additional genes relevant to other myeloid malignancies, refer to the [Myeloid Malignancies Mutation and Copy Number Variation Panel by Next Generation Sequencing](#) Test Fact Sheet.

For more information on ARUP's targeted rapid turnaround AML panel, which includes some of the genes on this panel, refer to the [Rapid Acute Myeloid Leukemia Targeted Therapy Mutation Panel](#) Test Fact Sheet.

For more information on ARUP's genomic microarray testing in oncology, refer to the [Cytogenomic Microarray - Oncology](#) Test Fact Sheet.

Limitations

Limit of Detection

- SNVs and small variants <24 base pairs (bp): 5% variant allele fraction (VAF)
- Variants >24 bp: May be detected at limit of detection (LOD) but analytic sensitivity may be reduced

Variants That Will Not Be Detected

- Variants in regions that are not included in the preferred transcript for the targeted genes
 - Refer to the [Genes Tested](#) table below for a full list of targeted regions and exclusions.
- Copy number variants (losses or gains)
- Loss of heterozygosity
- RNA variants
- Gene fusions, balanced translocations, and other structural variants
- Some variants may not be identified due to technical limitations in the presence of pseudogenes or in repetitive or homologous regions.

Variants That Will Not Be Reported

Benign or likely benign variants in the preferred transcript will not be reported.

Additional Limitations

- This panel is not intended to detect minimal residual disease (MRD).
- Interpretation of panel results may be impacted if the patient has had an undisclosed allogeneic bone marrow transplant or stem cell transplant.
- This panel does not distinguish between somatic and germline variants.

Analytic Sensitivity

Variant Class	Analytic Sensitivity (PPA) ^a Estimate (%)	Analytic Sensitivity (PPA) 95% Credibility Region ^a (%)
SNVs	>99	99.4-100.0
Indels, duplications, complex variants	98.5	96.3-99.5
<i>FLT3</i> -ITDs	96.6	89.4-99.3

^aGenes included on this panel are a subset of a larger methods-based validation from which the PPA values are derived.

ITDs, internal tandem duplications; PPA, positive percent agreement

Genes Tested

Genes Tested by Acute Myeloid Leukemia Mutation Panel by Next Generation Sequencing	
Gene	Preferred Transcript ^{a,b}
<i>ANKRD26</i>	NM_014915
<i>ASXL1</i>	NM_015338
<i>CEBPA</i>	NM_004364
<i>DDX41</i>	NM_016222
<i>DNMT3A</i>	NM_175629
<i>ETV6</i>	NM_001987

Gene	Preferred Transcript ^{a,b}
<i>FLT3</i>	NM_004119
<i>GATA2</i>	NM_032638
<i>IDH1</i>	NM_005896
<i>IDH2</i>	NM_002168
<i>KIT</i>	NM_000222
<i>KRAS</i>	NM_004985
<i>NPM1</i> ^c	NM_002520
<i>NRAS</i>	NM_002524
<i>RUNX1</i>	NM_001754
<i>TP53</i>	NM_000546
<i>WT1</i>	NM_024426

^aThis is the transcript number used for analyzing and reporting variants. The transcript version number may change periodically and thus is not listed here. The transcript with the version number will be included on the patient's report if a variant is detected in the gene.

^bNoncoding exons are not analyzed, except for regions containing known clinically relevant variants in the *ANKRD26* 5'UTR.

^cExon 1 is excluded due to technical limitations of the assay.

ARUP Laboratories is a nonprofit enterprise of the University of Utah and its Department of Pathology. 500 Chipeta Way, Salt Lake City, UT 84108
(800) 522-2787 | (801) 583-2787 | aruplab.com | arupconsult.com