

NEW TEST

[Click for Pricing](#)

T-Cell Clonality Screening, Tissue, by PCR

3021140, TCELLTISS

Specimen Requirements:

Patient Preparation:

Collect: Fresh tissue or formalin-fixed tumor tissue.

Specimen Preparation: Fresh Tissue: Freeze immediately. Transport 100 mg or 0.5-2.0 cm³ tissue
FFPE Tumor Tissue: Formalin-fixed (10 percent neutral buffered formalin) and paraffin-embedded tissue. Protect from excessive heat. Tissue block will be returned after testing. Transport tissue block or four 10-micron shavings in a tissue transport kit (ARUP Supply #47808) available online through eSupply using ARUP Connect(TM) or contact ARUP Client Services at 800-522-2787.

Transport Temperature: Fresh Tissue: Frozen on dry ice.
FFPE Tumor Tissue: Room temperature. Also acceptable: Refrigerated. Ship in cooled container during summer months.

Unacceptable Conditions: Tissue: FFPE specimens fixed in any fixative other than 10 percent neutral buffered formalin.

Remarks:

Stability: Fresh Tissue: Ambient: Unacceptable; Refrigerated: 2 hours; Frozen: 1 year
FFPE Tumor Tissue: Ambient: Indefinitely; Refrigerated: Indefinitely; Frozen: Unacceptable

Methodology: Qualitative Capillary Electrophoresis / Qualitative Polymerase Chain Reaction (PCR)

Note:

CPT Codes: 81342

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Refer to report.

Reference Interval:

HOTLINE NOTE: Refer to the Hotline Test Mix for interface build information.

