TEST CHANGE

Hepatitis A Virus Panel

0020597, HEP A PAN

Specimen Requirements:	
Patient Preparation:	Refer to individual components.
Collect:	Serum separator tube. Also acceptable: Lavender (EDTA) or pink (K2EDTA).
Specimen Preparation:	Separate serum or plasma from cells ASAP or within 2 hours of collection. Transfer 1 mL serum or plasma to an ARUP standard transport tube. (Min: 0.5 mL)
Transport Temperature:	Refrigerated.
Unacceptable Conditions:	Specimens containing particulate material. Severely hemolyzed or lipemic specimens. Heat-inactivated specimens. Specimens collected in citrate-based anticoagulant.
Remarks:	
Stability:	After separation from cells: Ambient: 48 hours~~Unacceptable~~; Refrigerated: 14 days~~1 week~~; Frozen: 30 days~~Indefinitely~~ (avoid repeated freeze/thaw cycles)
Methodology:	Qualitative Chemiluminescent Immunoassay (CLIA)
Note:	
CPT Codes:	86708; 86709
New York DOH Approval Status:	This test is New York DOH approved.
Interpretive Data:	

Reference Interval:

Test Number	Components	Reference Interval
	Hepatitis A Antibodies, Total	Negative
	Hepatitis A Antibody, IgM	Negative