NEW TEST

Click for Pricing

# Lipoprotein (a), Molarity

3021255, LP(A)

| | |
|---|---|
| Specimen Requirements: | |
| Patient Preparation: | |
| Collect: | Serum Separator Tube (SST) or Plasma Separator Tube (PST). Also acceptable: Green (Lithium Heparin), Lavender (EDTA), or Pink (K2EDTA). |
| Specimen Preparation: | Allow specimen to clot completely at room temperature. Separate from cells ASAP or within 2 hours of collection. Transfer 1 mL serum or plasma to an ARUP standard transport tube. (Min: 0.5 mL) |
| Transport Temperature: | Refrigerated. |
| Unacceptable Conditions: | Body fluids. |
| Remarks: | |
| Stability: | After separation from cells: Ambient: 24 hours; Refrigerated: 21 days; Frozen: 1 month. Freeze/Thaw: Up to 3 cycles. |
| Methodology: | Quantitative Immunoturbidimetry |
| Note: | |
| CPT Codes: | 83695 |
| New York DOH Approval Status: | This test is New York DOH approved. |

Interpretive Data:

ASCVD Relative Risk by Lp(a) Concentration

| Lp(a) (nmol/L) | Relative Risk of ASCVD |
|---|---|
| >430 | 4x |
| 350-429 | 3x |
| 250-349 | 2x |
| 125-249 | 1.4x |
| 75-124 | 1.2x |
| <75 | Reference |

Relative risks as shown in 2026 ACC/AHA Guideline on the Management of Dyslipidemia (Blumenthal et al., Circulation, 2026).

Reference Interval:

| Test Number | Components | Reference Interval |
|---|---|---|

| | | |
|---|---|---|
| | Lipoprotein (a) | <=74 |

HOTLINE NOTE: Refer to the Hotline Test Mix for interface build information.