

0025055 Heavy Metals Panel 6, Urine with Reflex to Arsenic Fractionated

HYMET 6

Reference Interval:

Test Number	Components	Reference Interval
0025000	Arsenic, Urine with Reflex to Fractionated	Effective November 13, 2017
		Test Number Components Reference Interval
		Arsenic, Urine- per volume 0.0- 34.9 µg/L (based on Biological Exposure Index)
		Arsenic, Urine- per24h 0.0- 49.9 µg/d
		Arsenic, Urine-ratio to CRT 0.0-29.9 ug/gCRT
0020734	Arsenic, Fractionated, Urine	Refer to report
	Creatinine, Urine - per 24h	Refer to report
0025040	Cadmium, Urine	Effective November 13, 2017
		Test Number Components Reference Interval
		Cadmium, Urine - per volume 0.0- 1.0 µg/L
		Cadmium, Urine - per 24h 0.0- 3.2 µg/d
		Cadmium, Urine - ratio to CRT 0.0- 3.2 µg/g CRT
Creatinine, Urine - per 24h Refer to report		
0020461	Copper, Urine	Effective November 13, 2017
		Test Number Components Reference Interval
		Copper, Urine- per volume 0.3-3.2 µg/dL
		Copper, Urine- per 24-h 3.0- 45.0 µg/d
		Copper, Urine-ratio to CRT 10.0-45.0 µg/g CRT
Creatinine, Urine - per 24h Refer to report		
0025060	Lead, Urine	Effective November 13, 2017
		Test Number Components Reference Interval
		Lead, Urine - per 24h 0.0- 8.1 µg/d
		Lead, Urine - per volume 0.0- 1.4 µg/L
		Lead Urine-ratio to CRT 0.0-1.4 ug/gCRT
Creatinine, Urine - per 24h Refer to report		
0025050	Mercury, Urine	Effective November 13, 2017
		Test Number Components Reference Interval
		Mercury, Urine - per 24h 0.0- 2.9 µg/d
		Mercury, Urine - per volume 0.0- 1.9 µg/L
		Mercury, Urine - ratio to CRT 0.0-20.0 µg/gCRT
Creatinine, Urine - per 24h Refer to report		
0020462	Zinc, Urine	Effective November 13, 2017
		Test Number Components Reference Interval
		Zinc, Urine-per volume 15.0-120.0 µg/dL
		Zinc, Urine-per 24h 150.0-1200.0 µg/d
		Zinc, Urine-ratio to CRT 110.0-750.0 µg/gCRT
Creatinine, Urine - per 24h Refer to report		

HOTLINE NOTE: There is a numeric map change associated with this test.

- Change the numeric map for component 0020101, Copper, Urine - per 24h from XXXXXX to XXXXXX.X.
- Change the numeric map for component 0025062, Lead, Urine - per volume from XXXXX to XXXXX.X.
- Change the numeric map for component 0025061, Lead, Urine - per 24h from XXXXX to XXXXX.X
- Change the numeric map for component 0025051, Mercury, Urine - per 24h from XXX to XXX.X.
- Change the numeric map for component 0020102, Zinc, Urine - per volume from XXXXXX to XXXXXX.X.
- Change the numeric map for component 0020103, Zinc, Urine - per 24h from XXXXXX to XXXXXX.X.