
0054441 Measles (Rubeola) Antibody, IgM, CSF

MEASLMCSF

Specimen Required: Collect: CSF.

Specimen Preparation: Transfer 0.5 mL CSF to an ARUP Standard Transport Tube. (Min: 0.2 mL)

Storage/Transport Temperature: Refrigerated. Also acceptable: Frozen.

Unacceptable Conditions: Specimen types other than CSF. Contaminated, heat-inactivated, or hemolyzed specimens.

Stability (collection to initiation of testing): Ambient: 8 hours; Refrigerated: 2 weeks; Frozen: 1 year