
0091161

Gamma-Hydroxybutyric Acid (GHB), Urine - Screen with Reflex to Confirmation/Quantitation

GHB U

Specimen Required: Collect: Random urine.

Specimen Preparation: Transfer 4 mL urine to an ARUP Standard Transport Tube. (Min: 2.4 mL)

Storage/Transport Temperature: Refrigerated. Also acceptable: Room temperature or frozen.

Stability (collection to initiation of testing): Ambient: 1 week; Refrigerated: 1 week; Frozen: 3 weeks