

---

**0098974      Angiotensin Converting Enzyme, CSF**

**ACE CSF**

**Specimen Required:** Collect: CSF.

Specimen Preparation: Separate from cells within 1 hour of collection. Transfer 1 mL CSF to an ARUP Standard Transport Tube.  
(Min: **0.5** mL)

Storage/Transport Temperature: Frozen.

Unacceptable Conditions: Hemolyzed or xanthochromic specimens.

Stability (collection to initiation of testing): Ambient: 4 hours; Refrigerated: 1 week; Frozen: 6 months