

Quarterly HOTLINE: Effective August 20, 2018

0099564 ***Strongyloides* Antibody, IgG by ELISA, Serum**

STRONGY

Specimen Required: Collect: Serum Separator Tube (SST) or Plain Red.

Specimen Preparation: Transfer 1 mL serum to an ARUP Standard Transport Tube. (Min. 0.3 mL)

Storage/Transport Temperature: Refrigerated.

Unacceptable Conditions: Bacterially contaminated, heat-inactivated, hemolyzed, icteric, or lipemic specimens.

Stability (collection to initiation of testing): After separation from cells: Ambient: 48 hours; Refrigerated: 2 weeks; Frozen: 1 year

Reference Interval:

Effective August 20, 2018

0.9 IV or less	Negative - No significant level of <i>Strongyloides</i> IgG antibody detected.
1.0 IV	Equivocal - The <i>Strongyloides</i> IgG antibody result is borderline and therefore inconclusive. Recommend retesting the patient in 2-4 weeks, if clinically indicated.
1.1 IV or greater	Positive - IgG antibodies to <i>Strongyloides</i> detected, which may suggest current or past infection.

HOTLINE NOTE: There is a numeric map change associated with this test.

Change the numeric map for component 0099564, *Strongyloides* Antibody, IgG By ELISA from XX.XX to XXX.X.