

Quarterly HOTLINE: Effective February 20, 2018

2001769

Bordetella pertussis Antibody, IgM by ELISA with Reflex to Immunoblot

BORDM

Reference Interval:

| Test Number | Components | Reference Interval | |
|-------------|---|--|--|
| | <i>Bordetella pertussis</i> Antibody, IgM by ELISA | Effective February 20, 2018 | |
| | | 0.9 IV or less | Negative - No significant level of detectable <i>B. pertussis</i> IgM antibody. |
| | | 1.0-1.1 IV | Equivocal - Repeat testing in 10-14 days may be helpful. |
| | | 1.2 IV or greater | Positive - IgM antibody to <i>B. pertussis</i> detected, which may indicate a current or recent exposure/immunization to <i>B. pertussis</i> . |
| 2004326 | <i>Bordetella pertussis</i> Antibody, IgM by Immunoblot | Effective February 19, 2013 | |
| | | Components | Reference Interval |
| | | <i>Bordetella pertussis</i> Ab, IgM by Immunoblot Interp | Negative |
| | | <i>B. pertussis</i> , IgM Immunoblot PT | Negative |
| | | <i>B. pertussis</i> , IgM Immunoblot FHA | Negative |

Note: If *Bordetella pertussis* Antibody, IgM by ELISA is 1.2 IV or greater, then *Bordetella pertussis* IgM Immunoblot testing will be added. Additional charges apply.

HOTLINE NOTE: There is a unit of measure change associated with this test. Change the unit of measure for component 2001783, *B. pertussis* Ab, IgM by ELISA from U/mL to IV.