

Quarterly HOTLINE: Effective February 20, 2018

2001774 *Bordetella pertussis* Antibodies, IgA and IgG by ELISA with Reflex to Immunoblot BORDPAN2

Reference Interval:

Test Number	Components	Reference Interval			
	<i>Bordetella pertussis</i> Antibody, IgA by ELISA with Reflex to Immunoblot	Test Number	Components	Reference Interval	
			<i>Bordetella pertussis</i> Antibody, IgA by ELISA	Effective February 20, 2018	
				0.9 IV or less	Negative - No significant level of detectable <i>B. pertussis</i> IgA antibody.
				1.0-1.1 IV	Equivocal - Repeat testing in 10-14 days may be helpful.
		1.2 IV or greater:	Positive - IgA antibody to <i>B. pertussis</i> detected, which may indicate a current or past exposure/immunization to <i>B. pertussis</i> .		
		2004316	<i>Bordetella pertussis</i> Antibody, IgA by Immunoblot	Negative	
2001768	<i>Bordetella pertussis</i> Antibody, IgG by ELISA with Reflex to Immunoblot	Test Number	Components	Reference Interval	
			<i>Bordetella pertussis</i> Antibody IgG by ELISA	Effective February 20, 2018	
				0.94 IV or less	Negative – No significant level of detectable <i>B. pertussis</i> IgG antibody.
				0.95-1.04 IV	Equivocal – Repeat testing in 10-14 days may be helpful.
		1.05 IV or greater	Positive – IgG antibody to <i>B. pertussis</i> detected, which may indicate a current or recent exposure/immunization to <i>B. pertussis</i> .		
		2004327	<i>Bordetella pertussis</i> Antibody, IgG by Immunoblot	Refer to report	

Note: If *Bordetella pertussis* Antibody, IgA by ELISA is 1.2 **IV** or greater, then *Bordetella pertussis* IgA Immunoblot testing will be added; if *Bordetella pertussis* Antibody, IgG by ELISA is 1.05 **IV** or greater, then *Bordetella pertussis* IgG Immunoblot testing will be added. Additional charges apply.

HOTLINE NOTE: There is a unit of measure and a numeric map change associated with this test.
 Change the unit of measure for component 2001770, *B. pertussis* Ab, IgA by ELISA from U/mL to **IV**.
 Change the unit of measure for component 2001782, *B. pertussis* Ab, IgG by ELISA from U/mL to **IV**.
 Change the numeric map for component 2001782, *B. pertussis* Ab, IgG by ELISA from X.X to **X.XX**.