

2005268 *Bordetella pertussis* Antibody, IgG by ELISA

BORDIGG

Reference Interval:

Effective February 20, 2018

0.94 IV or less	Negative – No significant level of detectable <i>B. pertussis</i> IgG antibody.
0.95-1.04 IV	Equivocal – Repeat testing in 10-14 days may be helpful.
1.05 IV or greater	Positive – IgG antibody to <i>B. pertussis</i> detected, which may indicate a current or recent exposure/immunization to <i>B. pertussis</i> .

Interpretive Data: See Compliance Statement D: www.aruplab.com/CS

HOTLINE NOTE: There is a unit of measure and a numeric map change associated with this test.
 Change the unit of measure for component 2001782, B. pertussis Ab, IgG by ELISA from U/mL to IV.
 Change the numeric map for component 2001782, B. pertussis Ab, IgG by ELISA from X.X to X.XX.