

HOTLINE: Effective August 19, 2019

---

**2005894**

**Allergen Panel, IgE by ImmunoCap ISAC**

**ISAC MICRO**

**Specimen Required:** Collect: Serum Separator Tube (SST). Also acceptable: Lavender (EDTA) or Green (Sodium or Lithium Heparin).

Specimen Preparation: Transfer 0.4 mL serum or plasma to an ARUP Standard Transport Tube. (Min: 0.25 mL)

Storage/Transport Temperature: Refrigerated. Also acceptable: Frozen.

Stability (collection to initiation of testing): Ambient: 24 hours; Refrigerated: 1 week; Frozen: 1 month

**CPT Code(s):** 86008 x112