
2008916

**Encephalitis Panel with Reflex to Herpes Simplex Virus Types 1 and 2
Glycoprotein G-Specific Antibodies, IgG, CSF**

ENCEPHCSF

Specimen Required: Collect: CSF.

Specimen Preparation: Transfer 3 mL CSF to an ARUP Standard Transport Tube. (Min: 1.05 mL)

Storage/Transport Temperature: Refrigerated.

Unacceptable Conditions: Serum or plasma. Contaminated, heat-inactivated, or hemolyzed specimens.

Stability (collection to initiation of testing): Ambient: 8 hours; Refrigerated: 2 weeks; Frozen: 1 year