

2010677

Parathyroid Hormone-Related Peptide (PTHrP) by LC-MS/MS, Plasma

PTHrP



Additional Technical Information



Specimen Collection and Handling

Specimen Required: Collect: Protease Inhibitor tube (PPACK; Phe-Pro-Arg-chlormethylketone) (ARUP supply #49662), available online through eSupply using ARUP Connect™ or contact ARUP Client Services at (800) 522-2787. **A winged collection set must be used.**
Specimen Preparation: Mix well. Separate from cells within 1 hour of collection. Transfer 1.5 mL plasma to an ARUP Standard Transport Tube. (Min: 0.7 mL)
Storage/Transport Temperature: Frozen. Separate specimens must be submitted when multiple tests are ordered.
Unacceptable Conditions: Grossly hemolyzed specimens.
Stability (collection to initiation of testing): After separation from cells: Ambient: Unacceptable; Refrigerated: Unacceptable; Frozen: 3 months