
2013715 LipoFit by NMR, Particle Count Only**NMRLIPFITP****Specimen Required:** Patient Prep: **Fasting specimen required.**

Collect: Greiner Bio-One Clot Activator Tube (ARUP supply #53483). Available online through eSupply using ARUP Connect™ or by contacting ARUP Client Services at (800) 522-2787. Also acceptable: Plain red.

Specimen Preparation: Gently invert tube to mix contents and allow to clot at room temperature. Separate serum from cells within 8 hours. Transfer 2 mL serum to an ARUP Standard Transport Tube. (Min: 1 mL)

Storage/Transport Temperature: Refrigerated.

Unacceptable Conditions: Plasma. Serum separator tubes other than Greiner Bio-One. **Non-fasting or lipemic specimens.**

Stability (collection to initiation of testing): Ambient: 2 days; Refrigerated: 1 month; Frozen: Unacceptable