

HOTLINE NOTE: There is a price change associated with this test. Please contact ARUP Client Services at (800) 522-2787 for additional information.

LipoFit by NMR Specimen Requirement Updates

ARUP Laboratories will be updating specimen requirements for LipoFit by NMR (2013716) and LipoFit by NMR, Particle Count Only (2013715) October 18, 2017.

Please note this change:

- Patient required to fast prior to specimen collection.

In order to reduce patient redraw due to analytical interference, please implement these updated specimen requirements as soon as possible.

Specimen collection reminder:

- Plain red-top tube or Greiner Bio-One Clot Activator tube (ARUP supply #53483) available online through eSupply using ARUP Connect™ or by contacting ARUP Client Services at (800) 522-2787.
 - Serum gel tubes, including BD (SST) tubes, are known to cause analytical interference.

More Information

Please refer inquiries regarding this notice to ARUP Client Services at (800) 522-2787, keyword: LipoFit

The change will be published with the February 20, 2018, quarterly Hotline notification.