

HOTLINE: Effective August 19, 2019

2013942 Zika Virus IgM Antibody Capture (MAC) by ELISA

ZIKA M

Interpretive Data: This assay is intended for in vitro diagnostic use under FDA Emergency Use Authorization (EUA) and has not been FDA cleared or approved. In compliance with this authorization, please visit <https://aruplab.com/zika> for more information and to access the applicable information sheets.

The possibility of false-positive or false-negative results must be considered. RT-PCR testing on both a serum and urine specimen is recommended by the Centers for Disease Control and Prevention (CDC) to rule out false-negative IgM results in patients experiencing symptoms for less than 2 weeks. Specimens collected for IgM testing greater than or equal to 2 weeks after symptom onset do not require any additional testing. For more information, please review the current clinical guidelines for Zika virus testing at: www.cdc.gov/zika/.

Note: If the result is “Presumptive Zika,” then Zika IgM Ab Capture (MAC) Confirmation (ARUP test code 3001904) will be added at no additional charge.

HOTLINE NOTE: There is a reflexive pattern change associated with this test.
Add reflex to 3001904, Zika IgM Ab Capture (MAC) Confirmation