
New Test **2014677** **Expanded Carrier Screen by Next Generation Sequencing with** **ECS SEQ FX**
Fragile X



Patient History for Expanded Carrier
Screening



Additional Technical Information

Methodology: Massively Parallel Sequencing/Polymerase Chain Reaction
Performed: Varies
Reported: Within 3 weeks

Specimen Required: Collect: Lavender (EDTA).
Specimen Preparation: Transport 4 mL whole blood. (Min: 1 mL)
Storage/Transport Temperature: Room temperature. Also acceptable: Refrigerated.
Remarks: Patient History form required.
Stability (collection to initiation of testing): Ambient: 1 week; Refrigerated: 1 week; Frozen: Unacceptable

Reference Interval: By report

CPT Code(s): 81479

New York DOH Approved.

HOTLINE NOTE: Refer to the Test Mix Addendum for interface build information.