

---

2014683

LeukoStrat CDx *FLT3* Mutation Detection by PCR

FLT3 CDX

**Specimen Required:** Collect: Green (Sodium Heparin).

Specimen Preparation: Transport 3 mL whole blood. (Min: 1 mL) **OR** 1 mL bone marrow. (Min: 0.25 mL) **Separate specimens must be submitted when multiple tests are ordered.**

Storage/Transport Temperature: Refrigerated.

Remarks: Specimen type required.

Unacceptable Conditions: **Specimens in Bone Marrow Transport Media.** Grossly hemolyzed or clotted specimens.

Stability (collection to initiation of testing): Ambient: 72 hours; Refrigerated: 1 week; Frozen: Unacceptable