

Quarterly HOTLINE: Effective November 5, 2018

New Test **3000394** **Malignancy Risk Assessment, Pelvic Mass, OVA1** **OVA1**

Methodology: Electrochemiluminescent Immunoassay (ECLIA)/Fixed-Rate-Time Nephelometry
Performed: Varies
Reported: 4-8 days

Specimen Required: Collect: Serum Separator Tube (SST).
Specimen Preparation: Transfer 2.2 mL serum to an ARUP Standard Transport Tube. (Min: 1.1 mL)
Storage/Transport Temperature: Refrigerated. Also acceptable: Frozen.
Stability (collection to initiation of testing): Ambient: Unacceptable; Refrigerated: 8 days; Frozen: 9 weeks

Reference Interval: By Report

Note: Biomarkers: CA-125 II, Apolipoprotein A1 (Apo A-1), Beta-2 Microglobulin (B2M), Transferrin, and Prealbumin.

CPT Code(s): 81503

New York DOH Approved.

HOTLINE NOTE: Refer to the Test Mix Addendum for interface build information.