

Client: Example Client ABC123
123 Test Drive
Salt Lake City, UT 84108
UNITED STATES

Physician: Doctor, Example

Patient: Patient, Example

DOB: 7/23/1982
Gender: Female
Patient Identifiers: 01234567890ABCD, 012345
Visit Number (FIN): 01234567890ABCD
Collection Date: 00/00/0000 00:00

Leukocyte Alkaline Phosphatase

ARUP test code 0049000

Leukocyte Alkaline Phosphatase

10 L (Ref Interval: 33-149)

TEST INFORMATION: Leukocyte Alkaline Phosphate (LAP) Stain Access complete set of age- and/ or gender specific reference intervals for this test in the ARUP Laboratory Test Directory (aruplab.com). The LAP stain determines cytoplasmic LAP enzyme activity in segmented neutrophils and bands. An increased/decreased LAP score can occur in various disease conditions. A frequent use of this test is to differentiate between Chronic Myelogenous Leukemia (CML) and a leukemoid reaction. CML has a decreased LAP score; leukemoid reaction has an increased LAP score.

Methodology: Enzyme cytochemistry The LAP enzyme present in neutrophils hydrolyzes the substrate Naphthol AS-BI Phosphate. The hydrolyzed product couples with the dye (fast red violet salt L.B.) to form an insoluble red-brown precipitate in the cytoplasm. Smears are counterstained with hematoxylin. The smears are examined microscopically. Results of this test must be interpreted within the clinical context. This test should not be used alone to determine diagnosis. Test developed by ARUP Laboratories.

VERIFIED/REPORTED DATES

Procedure	Accession	Collected	Received	Verified/Reported
Leukocyte Alkaline Phosphatase	21-284-115334	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00

END OF CHART

H=High, L=Low, *=Abnormal, C=Critical

Unless otherwise indicated, testing performed at: