

Client: Example Client ABC123
123 Test Drive
Salt Lake City, UT 84108
UNITED STATES

Physician: Doctor, Example

Patient: Patient, Example

DOB: 11/2/1990
Gender: Female
Patient Identifiers: 01234567890ABCD, 012345
Visit Number (FIN): 01234567890ABCD
Collection Date: 00/00/0000 00:00

Chlamydia Antibody Panel, IgG & IgM by IFA

ARUP test code 0065100

C. pneumoniae IgM Titer <1:20 (Ref Interval: <1:20)

C. trachomatis IgM Titer <1:20 (Ref Interval: <1:20)

C. psittaci IgM Titer <1:20 (Ref Interval: <1:20)

C. pneumoniae IgG Titer **1:512** * (Ref Interval: <1:64)

C. trachomatis IgG Titer **1:512** * (Ref Interval: <1:64)

C. psittaci IgG Titer **1:64** * (Ref Interval: <1:64)

INTERPRETIVE INFORMATION: C. psittaci IgG Titer

The chlamydia antibody test contains both species- and genus-specific antigens, and serological cross-reactions may be seen in both acute and convalescent samples (less than 1:128). A C. pneumoniae-specific reaction will exhibit titers twofold or greater than titers observed with the C. trachomatis or C. psittaci serology. Any IgG titer may indicate past exposure to that particular species. IgG titers in recently infected individuals are typically greater than or equal 1:512.

The chlamydia microimmunofluorescent assay slides utilize C. psittaci, C. pneumoniae, and nine serotypes of C. trachomatis. The LGV strains of C. trachomatis are not included in this assay.

Test developed and characteristics determined by ARUP Laboratories. See Compliance Statement A: aruplab.com/CS

H=High, L=Low, *=Abnormal, C=Critical

VERIFIED/REPORTED DATES				
Procedure	Accession	Collected	Received	Verified/Reported
C. pneumoniae IgM Titer	20-289-113519	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
C. trachomatis IgM Titer	20-289-113519	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
C. psittaci IgM Titer	20-289-113519	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
C. pneumoniae IgG Titer	20-289-113519	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
C. trachomatis IgG Titer	20-289-113519	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
C. psittaci IgG Titer	20-289-113519	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00

END OF CHART

H=High, L=Low, *=Abnormal, C=Critical

Unless otherwise indicated, testing performed at:

ARUP LABORATORIES | 800-522-2787 | aruplab.com
500 Chipeta Way, Salt Lake City, UT 84108-1221
Tracy I. George, MD, Laboratory Director

Patient: Patient, Example
ARUP Accession: 20-289-113519
Patient Identifiers: 01234567890ABCD, 012345
Visit Number (FIN): 01234567890ABCD
Page 2 of 2 | Printed: 10/20/2020 10:54:03 AM
4848