

Client: Example Client ABC123  
123 Test Drive  
Salt Lake City, UT 84108  
UNITED STATES

Physician: Doctor, Example

**Patient: Patient, Example**

**DOB:** 8/14/1985  
**Gender:** Female  
**Patient Identifiers:** 01234567890ABCD, 012345  
**Visit Number (FIN):** 01234567890ABCD  
**Collection Date:** 01/01/2017 12:34

**Drug Panel 7, Urine - Screen with Reflex to Confirmation/Quantitation**

ARUP test code 0092184

Amphetamines, Urn, Screen	Negative ng/mL	(Ref Interval: Cutoff 300)
Barbiturates, Urn, Screen	Negative ng/mL	(Ref Interval: Cutoff 200)
Benzodiazepines, Urn, Screen	Negative ng/mL	(Ref Interval: Cutoff 200)
Cocaine, Urn, Screen	Negative ng/mL	(Ref Interval: Cutoff 150)
Creatinine, Urine	29.2 mg/dL	(Ref Interval: 20.0-400.0)
Opiates, Urn, Screen	Negative ng/mL	(Ref Interval: Cutoff 300)
Phencyclidine, Urn, Screen	Negative ng/mL	(Ref Interval: Cutoff 25)
THC, Urn, Screen	Negative ng/mL	(Ref Interval: Cutoff 20)
CDASU 7 Comments	See Note	

**H=High, L=Low, \*=Abnormal, C=Critical**

Unless otherwise indicated, testing performed at:

INTERPRETIVE INFORMATION: Drug Panel 7, Urn, Scrn w/Rflx to Conf

The absence of expected drug(s) and/or drug metabolite(s) may indicate non-compliance, inappropriate timing of specimen collection relative to drug administration, poor drug absorption, diluted/adulterated urine, or limitations of testing. The concentration at which the screening test can detect a drug or metabolite varies within a drug class. Specimens for which drugs or drug classes are detected by the screen are reflexed to a second, more specific technology (GC/MS and/or LC-MS/MS). The concentration value must be greater than or equal to the cutoff to be reported as positive. Interpretive questions should be directed to the laboratory.

For medical purposes only; not valid for forensic use.

Oxycodone results are reported with the opiates results. MDMA results are reported with the amphetamines results.

VERIFIED/REPORTED DATES

Procedure	Accession	Collected	Received	Verified/Reported
Amphetamines, Urn, Screen	20-067-108064	3/7/2020 10:54:00 PM	3/9/2020 11:44:42 AM	3/9/2020 6:38:00 PM
Barbiturates, Urn, Screen	20-067-108064	3/7/2020 10:54:00 PM	3/9/2020 11:44:42 AM	3/9/2020 6:38:00 PM
Benzodiazepines, Urn, Screen	20-067-108064	3/7/2020 10:54:00 PM	3/9/2020 11:44:42 AM	3/9/2020 6:38:00 PM
Cocaine, Urn, Screen	20-067-108064	3/7/2020 10:54:00 PM	3/9/2020 11:44:42 AM	3/9/2020 6:38:00 PM
Creatinine, Urine	20-067-108064	3/7/2020 10:54:00 PM	3/9/2020 11:44:42 AM	3/9/2020 6:38:00 PM
Opiates, Urn, Screen	20-067-108064	3/7/2020 10:54:00 PM	3/9/2020 11:44:42 AM	3/9/2020 6:38:00 PM
Phencyclidine, Urn, Screen	20-067-108064	3/7/2020 10:54:00 PM	3/9/2020 11:44:42 AM	3/9/2020 6:38:00 PM
THC, Urn, Screen	20-067-108064	3/7/2020 10:54:00 PM	3/9/2020 11:44:42 AM	3/9/2020 6:38:00 PM
CDASU 7 Comments	20-067-108064	3/7/2020 10:54:00 PM	3/9/2020 11:44:42 AM	3/9/2020 6:38:00 PM

END OF CHART

H=High, L=Low, \*=Abnormal, C=Critical

Unless otherwise indicated, testing performed at: