

Client: Example Client ABC123
123 Test Drive
Salt Lake City, UT 84108
UNITED STATES

Physician: Doctor, Example

Patient: Patient, Example

DOB: 12/23/1958
Gender: Male
Patient Identifiers: 01234567890ABCD, 012345
Visit Number (FIN): 01234567890ABCD
Collection Date: 00/00/0000 00:00

Drug Panel 9, Urine - Screen with Reflex to Confirmation/Quantitation

ARUP test code 0092186

Amphetamines, Urn, Screen	Negative ng/mL	(Ref Interval: Cutoff 300)
Barbiturates, Urn, Screen	Negative ng/mL	(Ref Interval: Cutoff 200)
Benzodiazepines, Urn, Screen	Negative ng/mL	(Ref Interval: Cutoff 200)
Cocaine, Urn, Screen	Negative ng/mL	(Ref Interval: Cutoff 150)
Creatinine, Urine	95.5 mg/dL	(Ref Interval: 20.0-400.0)
Methadone, Urn, Screen	Negative ng/mL	(Ref Interval: Cutoff 150)
Opiates, Urn, Screen	Negative ng/mL	(Ref Interval: Cutoff 300)
Phencyclidine, Urn, Screen	Negative ng/mL	(Ref Interval: Cutoff 25)
Propoxyphene, Urn, Screen	Negative ng/mL	(Ref Interval: Cutoff 300)
THC, Urn, Screen	Negative ng/mL	(Ref Interval: Cutoff 50)
CDASU 9 Comments	See Note	

H=High, L=Low, *=Abnormal, C=Critical

Unless otherwise indicated, testing performed at:

INTERPRETIVE INFORMATION: Drug Panel 9, Urn, Scrn w/Rflx to Conf

The absence of expected drug(s) and/or drug metabolite(s) may indicate non-compliance, inappropriate timing of specimen collection relative to drug administration, poor drug absorption, diluted/adulterated urine, or limitations of testing. The concentration at which the screening test can detect a drug or metabolite varies within a drug class. Specimens for which drugs or drug classes are detected by the screen are reflexed to a second, more specific technology (GC/MS and/or LC-MS/MS). The concentration value must be greater than or equal to the cutoff to be reported as positive. Interpretive questions should be directed to the laboratory.

Oxycodone results are reported with the opiates results. MDMA results are reported with the amphetamines results. The following opioids are not detected in this test: fentanyl, buprenorphine, meperidine, tramadol, and tapentadol. A comprehensive panel that includes these opioids is available or individual opioid testing can be ordered. Refer to aruplab.com for test information.

For medical purposes only; not valid for forensic use.

VERIFIED/REPORTED DATES

Procedure	Accession	Collected	Received	Verified/Reported
Amphetamines, Urn, Screen	24-081-401173	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
Barbiturates, Urn, Screen	24-081-401173	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
Benzodiazepines, Urn, Screen	24-081-401173	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
Cocaine, Urn, Screen	24-081-401173	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
Creatinine, Urine	24-081-401173	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
Methadone, Urn, Screen	24-081-401173	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
Opiates, Urn, Screen	24-081-401173	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
Phencyclidine, Urn, Screen	24-081-401173	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
Propoxyphene, Urn, Screen	24-081-401173	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
THC, Urn, Screen	24-081-401173	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
CDASU 9 Comments	24-081-401173	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00

END OF CHART

H=High, L=Low, *=Abnormal, C=Critical

Unless otherwise indicated, testing performed at: