

Client: Example Client ABC123  
123 Test Drive  
Salt Lake City, UT 84108  
UNITED STATES

Physician: Doctor, Example

**Patient: Patient, Example**

**DOB:** 3/26/1983  
**Gender:** Female  
**Patient Identifiers:** 01234567890ABCD, 012345  
**Visit Number (FIN):** 01234567890ABCD  
**Collection Date:** 00/00/0000 00:00

**Bordetella pertussis Antibody, IgM by ELISA with Reflex to Immunoblot**

ARUP test code 2001769

B. pertussis Ab, IgM by ELISA

0.9 IV (Ref Interval: <=1.1)

When Bordetella pertussis Antibodies by ELISA testing is negative, no further testing is performed.

INTERPRETIVE INFORMATION: B. pertussis Ab, IgM w/Reflex

- 0.9 IV or less: Negative - No significant level of detectable B. pertussis IgM antibody.
- 1.0-1.1 IV: Equivocal - Repeat testing in 10-14 days may be helpful.
- 1.2 IV or greater: Positive - IgM antibody to B. pertussis detected, which may indicate a current or recent exposure/immunization to B. pertussis.

Recommend that treatment decisions be based on the result of the B. pertussis IgM immunoblot test instead of the ELISA test. B. pertussis IgM test by ELISA may produce false-positive results.

This test was developed and its performance characteristics determined by ARUP Laboratories. It has not been cleared or approved by the US Food and Drug Administration. This test was performed in a CLIA certified laboratory and is intended for clinical purposes.

VERIFIED/REPORTED DATES

Procedure	Accession	Collected	Received	Verified/Reported
B. pertussis Ab, IgM by ELISA	23-017-138337	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00

END OF CHART

H=High, L=Low, \*=Abnormal, C=Critical

Unless otherwise indicated, testing performed at: