

Client: Example Client ABC123
123 Test Drive
Salt Lake City, UT 84108
UNITED STATES

Physician: Doctor, Example

Patient: Patient, Example

DOB: 6/15/1949
Gender: Male
Patient Identifiers: 01234567890ABCD, 012345
Visit Number (FIN): 01234567890ABCD
Collection Date: 00/00/0000 00:00

Inhibitor Assay, PT with Reflex to PT 1:1 Mix

ARUP test code 2003260

PT, Inhibitor Screen **24.2 sec H** (Ref Interval: 12.0-15.5)

PT 1:1 mixing study added.
A PT-based inhibitor is detected.
The PT is prolonged and does not correct in a 1:1 mix with pooled normal plasma. Possibilities include a specific factor inhibitor, high heparin concentration, direct thrombin inhibitor, factor Xa inhibitor, dysfibrinogenemia, or monoclonal proteins. This pattern can occasionally be seen in the presence of multiple or severe factor deficiencies or with rare lupus anticoagulants. Please correlate clinically.

PT, Inhibitor Screen, 1:1 Mix

ARUP test code 2003264

PT, Inhibitor Screen, 1:1 Mix **18.8 sec H** (Ref Interval: 12.0-15.5)

VERIFIED/REPORTED DATES

Procedure	Accession	Collected	Received	Verified/Reported
PT, Inhibitor Screen	23-037-402056	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
PT, Inhibitor Screen, 1:1 Mix	23-037-402056	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00

END OF CHART

H=High, L=Low, *=Abnormal, C=Critical