

Client: Example Client ABC123
123 Test Drive
Salt Lake City, UT 84108
UNITED STATES

Physician: Doctor, Example

Patient: Patient, Example

DOB: 1/7/1945
Gender: Male
Patient Identifiers: 01234567890ABCD, 012345
Visit Number (FIN): 01234567890ABCD
Collection Date: 00/00/0000 00:00

Inhibitor Assay, PTT with Reflex to PTT 1:1 Mix, with Reflex to 1-Hour Incubation

ARUP test code 2003266

PTT, Inhibitor Screen **55 sec H** (Ref Interval: 32-48)

PTT 1:1 mixing study added.

A PTT-based inhibitor is not detected. The PTT is prolonged, corrects in a 1:1 mix with pooled normal plasma, and does not demonstrate time-dependence in the incubated mixing study. This pattern is most commonly due to coagulation factor deficiency.

Partial Thromboplastin Time 1:1 Mix (performed if PTT > 48 seconds)

ARUP test code 2003270

PTT, Inhibitor Screen, 1:1 Mix **43 sec** (Ref Interval: 32-48)

Incubated mixing study added.

Partial Thromboplastin Time 1 hour

ARUP test code 2003272

PTT, Inhibitor Screen, 1-Hour **58 sec**

PTT, Inhibitor Screen, 1-Hour, 1:1 Mix **43 sec**

H=High, L=Low, *=Abnormal, C=Critical

VERIFIED/REPORTED DATES				
Procedure	Accession	Collected	Received	Verified/Reported
PTT, Inhibitor Screen	23-030-134945	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
PTT, Inhibitor Screen, 1:1 Mix	23-030-134945	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
PTT, Inhibitor Screen, 1-Hour	23-030-134945	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
PTT, Inhibitor Screen, 1-Hour, 1:1 Mix	23-030-134945	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00

END OF CHART

H=High, L=Low, *=Abnormal, C=Critical

Unless otherwise indicated, testing performed at:

ARUP LABORATORIES | 800-522-2787 | aruplab.com
500 Chipeta Way, Salt Lake City, UT 84108-1221
Jonathan R. Genzen, MD, PhD, Laboratory Director

Patient: Patient, Example
ARUP Accession: 23-030-134945
Patient Identifiers: 01234567890ABCD, 012345
Visit Number (FIN): 01234567890ABCD
Page 2 of 2 | Printed: 3/7/2023 12:18:08 PM
4848