

Client: Example Client ABC123  
123 Test Drive  
Salt Lake City, UT 84108  
UNITED STATES

Physician: Doctor, Example

**Patient: Patient, Example**

**DOB:** 11/15/1992  
**Sex:** Female  
**Patient Identifiers:** 01234567890ABCD, 012345  
**Visit Number (FIN):** 01234567890ABCD  
**Collection Date:** 01/01/2017 12:34

**Ethyl Glucuronide Screen with Reflex to Confirmation, Urine**

ARUP test code 2007912

Ethyl Glucuronide Screen w/Rflx, Urine

Negative ng/mL (Ref Interval: Cutoff 500)

INTERPRETIVE INFORMATION: Ethyl Glucuronide Screen with Reflex to Confirmation, Urine

Ethyl glucuronide is a direct metabolite of ethanol and can be detected up to 80 hours in urine after ethanol ingestion. The cutoff for positive by immunoassay is set at 500 ng/mL. A positive result will be confirmed by liquid chromatography tandem mass spectrometry (LC-MS/MS).

VERIFIED/REPORTED DATES

Procedure	Accession	Collected	Received	Verified/Reported
Ethyl Glucuronide Screen w/Rflx, Urine	22-081-103868	3/22/2022 11:35:00 AM	3/22/2022 4:09:28 PM	3/22/2022 8:14:00 PM

END OF CHART

H=High, L=Low, \*=Abnormal, C=Critical

Unless otherwise indicated, testing performed at:

ARUP LABORATORIES | 800-522-2787 | aruplab.com  
500 Chipeta Way, Salt Lake City, UT 84108-1221  
Jonathan R. Genzen, MD, PhD, Laboratory Director

Patient: Patient, Example  
ARUP Accession: 22-081-103868  
Patient Identifiers: 01234567890ABCD, 012345  
Visit Number (FIN): 01234567890ABCD  
Page 1 of 1 | Printed: 9/13/2022 7:23:33 AM