

Client: Example Client ABC123 123 Test Drive Salt Lake City, UT 84108

UNITED STATES

Physician: Doctor, Example

**Patient: Patient, Example** 

**DOB** 2/9/1986

**Gender:** Male

**Patient Identifiers:** 01234567890ABCD, 012345

**Visit Number (FIN):** 01234567890ABCD **Collection Date:** 00/00/0000 00:00

## Fentanyl, Urine Screen with Reflex to Quantitation

ARUP test code 2012284

Fentanyl, Urn, Screen

PresumptivePOS ng/mL

(Ref Interval: Cutoff 1)

If the screen is PresumptivePos, then confirmation testing by mass spectrometry will be added. Additional charges will apply.

Screen, Urine Interpretation

See Note

INTERPRETIVE INFORMATION:

The absence of expected drug(s) and/or drug metabolite(s) may indicate non-compliance, inappropriate timing of specimen collection relative to drug administration, poor drug absorption, diluted/adulterated urine, or limitations of testing. The concentration at which the screening test can detect a drug or metabolite varies. Specimens for which drugs or drug classes are detected by the screen are reflexed to a second, more specific technology (GC/MS and/or LC-MS/MS). The concentration value must be greater than or equal to the cutoff to be reported as positive. Interpretive questions should be directed to the laboratory.

For medical purposes only; not valid for forensic use.

## Fentanyl and Metabolite, Urine, Quantitative

ARUP test code 0092570

Fentanyl, Urn, Quant

<1.0 ng/mL

H=High, L=Low, \*=Abnormal, C=Critical

4848



INTERPRETIVE INFORMATION: Fentanyl and Metabolite, Urine

Methodology: Quantitative Liquid Chromatography-Tandem Mass Spectrometry

Positive cutoff: 1.0 ng/mL

For medical purposes only; not valid for forensic use.

The absence of expected drug(s) and/or drug metabolite(s) may indicate non-compliance, inappropriate timing of specimen collection relative to drug administration, poor drug absorption, diluted/adulterated urine, or limitations of testing. The concentration value must be greater than or equal to the cutoff to be reported as positive. Interpretive questions should be directed to the laboratory.

This test was developed and its performance characteristics determined by ARUP Laboratories. It has not been cleared or approved by the US Food and Drug Administration. This test was performed in a CLIA certified laboratory and is intended for clinical purposes.

Norfentanyl, Urn, Quant

19.3 ng/mL

| VERIFIED/REPORTED DATES      |               |                  |                  |                   |
|------------------------------|---------------|------------------|------------------|-------------------|
| Procedure                    | Accession     | Collected        | Received         | Verified/Reported |
| Fentanyl, Urn, Screen        | 24-081-119680 | 00/00/0000 00:00 | 00/00/0000 00:00 | 00/00/0000 00:00  |
| Screen, Urine Interpretation | 24-081-119680 | 00/00/0000 00:00 | 00/00/0000 00:00 | 00/00/0000 00:00  |
| Fentanyl, Urn, Quant         | 24-081-119680 | 00/00/0000 00:00 | 00/00/0000 00:00 | 00/00/0000 00:00  |
| Norfentanyl, Urn, Quant      | 24-081-119680 | 00/00/0000 00:00 | 00/00/0000 00:00 | 00/00/0000 00:00  |

END OF CHART

H=High, L=Low, \*=Abnormal, C=Critical

Patient: Patient, Example ARUP Accession: 24-081-119680 Patient Identifiers: 01234567890ABCD, 012345 Visit Number (FIN): 01234567890ABCD Page 2 of 2 | Printed: 4/1/2024 3:58:56 PM

4848