

Client: Example Client ABC123 123 Test Drive

Salt Lake City, UT 84108 UNITED STATES

Physician: Doctor, Example

Patient: Patient, Example

DOB 4/3/1959 **Gender:** Male

Patient Identifiers: 01234567890ABCD, 012345

Visit Number (FIN): 01234567890ABCD **Collection Date:** 00/00/0000 00:00

Zolpidem Urine with Reflex to Quantitation

ARUP test code 2012300

Zolpidem, Urn, Screen

PresumptivePOS ng/mL

(Ref Interval: Cutoff 20)

If the screen is PresumptivePos, then confirmation testing by mass spectrometry will be added. Additional charges will apply.

Screen, Urine Interpretation

See Note

INTERPRETIVE INFORMATION:

The absence of expected drug(s) and/or drug metabolite(s) may indicate non-compliance, inappropriate timing of specimen collection relative to drug administration, poor drug absorption, diluted/adulterated urine, or limitations of testing. The concentration at which the screening test can detect a drug or metabolite varies. Specimens for which drugs or drug classes are detected by the screen are reflexed to a second, more specific technology (GC/MS and/or LC-MS/MS). The concentration value must be greater than or equal to the cutoff to be reported as positive. Interpretive questions should be directed to the laboratory.

For medical purposes only; not valid for forensic use.

Zolpidem, Urine, Quantitative

ARUP test code 2012319

Zolpidem, Urn, Quant

78 ng/mL

Consistent with use of a drug containing zolpidem, such as Ambien.

H=High, L=Low, *=Abnormal, C=Critical

4848



INTERPRETIVE INFORMATION: Zolpidem, Urine, Quantitative

Methodology: Quantitative Liquid Chromatography-Tandem Mass

Spectrometry
Drugs covered: zolpidem.

Positive cutoff: 20 ng/mL

For medical purposes only; not valid for forensic use.

The absence of expected drug may indicate non-compliance, inappropriate timing of specimen collection relative to drug administration, poor drug absorption, diluted/adulterated urine, or limitations of testing. The concentration value must be greater than or equal to the cutoff to be reported as positive. Interpretive questions should be directed to the laboratory.

This test was developed and its performance characteristics determined by ARUP Laboratories. It has not been cleared or approved by the US Food and Drug Administration. This test was performed in a CLIA certified laboratory and is intended for clinical purposes.

VERIFIED/REPORTED DATES				
Procedure	Accession	Collected	Received	Verified/Reported
Zolpidem, Urn, Screen	23-341-127004	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
Screen, Urine Interpretation	23-341-127004	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
Zolpidem, Urn, Quant	23-341-127004	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00

END OF CHART

H=High, L=Low, *=Abnormal, C=Critical