

Client: Example Client ABC123 123 Test Drive Salt Lake City, UT 84108 UNITED STATES

Physician: Doctor, Example

Patient: Patient, Example

DOB 4/19/1974

Gender: Male

Patient Identifiers: 01234567890ABCD, 012345

Visit Number (FIN): 01234567890ABCD **Collection Date:** 00/00/0000 00:00

Coccidioides Antibodies, IgM and IgG by Immunoassay, CSF

ARUP test code 3000057

Coccidioides Antibody IgG ELISA, CSF

0.0 IV (Ref Interval: <=0.2)

INTERPRETIVE INFORMATION: Coccidioides Antibody IgG ELISA,

0.2 IV or less:

Negative - No significant level of Coccidioides IgG antibody detected.

0.3 IV or greater:

Positive - Presence of IgG antibody to Coccidioides detected, suggestive of current or past infection.

IgG antibodies usually appear by the third week of infection and

may persist for years. Both tube precipitin (TP) and CF antigens are represented in the ELISA tests.

Coccidioides Antibody IgM ELISA, CSF

0.0 IV (Ref Interval: <=0.2)

INTERPRETIVE INFORMATION: Coccidioides Antibody Igm ELISA,

0.2 IV or less:

Negative - No significant level of

Coccidioides IgM antibody detected.

0.3 IV or greater:

Positive - Presence of IgM antibody to Coccidioides detected, suggestive

of current or past infection.

In most symptomatic patients, IgM antibodies usually appear by the second week of infection and disappear by the fourth month. Both tube precipitin (TP) and CF antigens are represented in the

ELISA tests.

H=High, L=Low, *=Abnormal, C=Critical



VERIFIED/REPORTED DATES				
Procedure	Accession	Collected	Received	Verified/Reported
Coccidioides Antibody IgG ELISA, CSF	23-039-153681	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
Coccidioides Antibody IgM ELISA, CSF	23-039-153681	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00

END OF CHART

H=High, L=Low, *=Abnormal, C=Critical