

Client: Example Client ABC123
123 Test Drive
Salt Lake City, UT 84108
UNITED STATES

Physician: Doctor, Example

Patient: Patient, Example

DOB: Unknown
Gender: Male
Patient Identifiers: 01234567890ABCD, 012345
Visit Number (FIN): 01234567890ABCD
Collection Date: 00/00/0000 00:00

RhE/e (RHCE) Antigen Genotyping

ARUP test code 3002003

RHE GENO Specimen

Amniotic fluid

RhEe Genotype

E/e

Indication for testing: Determine fetal RhEe genotype to assess risk for alloimmune hemolytic disease of the fetus and newborn (HDFN).

Fetal RhEe genotype: E/e

Interpretation: One copy of the RHCE*3 (E) allele and one copy of the RHCE*5 (e) allele were detected in this prenatal sample. This genotype is predictive of an RhE+e+ phenotype in this fetus.

H=High, L=Low, *=Abnormal, C=Critical

BACKGROUND INFORMATION: RhE/e (RHCE) Antigen Genotyping

CHARACTERISTICS: Erythrocyte alloimmunization may result in hemolytic transfusion reactions or hemolytic disease of the fetus and newborn (HDFN).
E ANTIGEN FREQUENCY: 0.29 Caucasians, 0.22 African Americans, 0.39 Asians.
e ANTIGEN FREQUENCY: 0.98 Caucasians, 0.98 African Americans, 0.96 Asians.
INHERITANCE: Co-dominant.
CAUSE: Antigen-antibody mediated red-cell hemolysis between donor/recipient or transferred maternal antibodies.
POLYMORPHISM TESTED: Rh blood group RHCE*3 (E), RHCE*5 (e): c.676G>C; p.Ala226Pro.
CLINICAL SENSITIVITY: 99 percent.
METHODOLOGY: Immucor PreciseType(TM) HEA Molecular BeadChip which is FDA-approved for clinical testing.
ANALYTIC SENSITIVITY AND SPECIFICITY: 99 percent.
LIMITATIONS: Bloody amniotic fluid samples may give false-negative results because of maternal cell contamination. Rare nucleotide changes leading to altered or partial antigen expression and null phenotypes are not detected by this assay. This assay is occasionally limited in predicting genotype due to extreme variation in the Rh locus. False-negative Rho predictions may result due to RHCE-D-CE fusion genes. Patients who have had hematopoietic stem cell transplants may have inconclusive results on this test. Abnormal signal intensities may result in indeterminate genotyping results.

For quality assurance purposes, ARUP Laboratories will confirm the above result at no charge following delivery. Order Confirmation of Fetal Testing and include a copy of the original fetal report (or the mother's name and date of birth) with the test submission. Please contact an ARUP genetic counselor at (800) 242-2787 extension 2141 prior to specimen submission.

VERIFIED/REPORTED DATES

Procedure	Accession	Collected	Received	Verified/Reported
RHE GENO Specimen	19-326-106032	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
RhEe Genotype	19-326-106032	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00

END OF CHART

H=High, L=Low, *=Abnormal, C=Critical