

Client: Example Client ABC123
123 Test Drive
Salt Lake City, UT 84108
UNITED STATES

Physician: Doctor, Example

Patient: Patient, Example

DOB: 3/13/1970
Gender: Female
Patient Identifiers: 01234567890ABCD, 012345
Visit Number (FIN): 01234567890ABCD
Collection Date: 00/00/0000 00:00

B Cell Subset Analysis

ARUP test code 3002216

CD19+ B cells %	7.3 % Lymphs	(Ref Interval: 6.4-22.0)
CD19+ B cells	134 cells/uL	(Ref Interval: 110-450)
CD20+ %	99.6 % of CD19	(Ref Interval: 96.0-100.0)
CD20+	133 cells/uL	(Ref Interval: 110-450)
Total Memory CD27+ %	10.6 % of CD19	(Ref Interval: 10.0-33.0)
Total Memory CD27+	14 cells/uL L	(Ref Interval: 23-110)
Non switched CD27+IgD+IgM+ %	4.1 % of CD19	(Ref Interval: 2.4-15.0)
Non switched CD27+IgD+IgM+	5 cells/uL	(Ref Interval: 5-46)
Class-switched CD27+IgD-IgM- %	5.1 % of CD19	(Ref Interval: 5.1-22.0)
Class-switched CD27+IgD-IgM-	7 cells/uL L	(Ref Interval: 11-61)
Transitional CD38+IgM+ %	5.0 % of CD19	(Ref Interval: 0.7-5.9)
Transitional CD38+IgM+	7 cells/uL	(Ref Interval: 1-17)

H=High, L=Low, *=Abnormal, C=Critical

Unless otherwise indicated, testing performed at:

Plasmablasts CD38+IgM- %	0.6 % of CD19	(Ref Interval: 0.4-4.1)
--------------------------	---------------	-------------------------

Plasmablasts CD38+IgM-	1 cells/uL	(Ref Interval: 1-8)
------------------------	------------	---------------------

Activated CD21low CD38- %	2.0 % of CD19	(Ref Interval: 1.2-9.0)
---------------------------	---------------	-------------------------

Activated CD21low CD38-	3 cells/uL	(Ref Interval: 3-26)
-------------------------	------------	----------------------

INTERPRETIVE INFORMATION: B Cell Subset Analysis

This panel identifies B-cell dysregulation. B-cells start development in the bone marrow (stem-cell, pro-B, pre-B), then transition to the spleen and lymph nodes where some mature by acquiring CD27 and switching immunoglobulin class from IgD and IgM to IgG or IgA. Class-switched B-cells may further progress to plasmablasts and finally plasma cells. Different disorders may block different parts of this pathway, disrupting immunoglobulin production.

This panel can also be used to monitor B-cell reconstitution after bone marrow transplantation or targeted B-cell depletion therapy.

This panel can assist in the diagnosis and subclassification of Common Variable Immune Deficiency (CVID). CVID is a heterogeneous group of disorders characterized by low antibody production, defective antibody responses, and recurrent infections. Most cases of CVID have a severe reduction in class switched memory B-cells (CD27+, IgD-, IgM-) that correlates with granulomatous disease. Many also have an expanded population of CD21low, CD38low B-cells that correlates with splenomegaly. Increased transitional B-cells (CD38+, IgM+) in CVID correlates with lymphadenopathy. Most CVID patients have a low percentage of plasmablasts (CD38+, IgM-) that has a correlation with autoimmune cytopenia.

Class switched memory B-cells are also low in ALPS, but are typically increased in SLE and infection.

Please note: Reference intervals for CD20+ B-cells were not established for patients less than 16 years of age. For all other B-cell subsets, reference intervals for populations younger than 16-years are adopted from literature. Piatosa B, Wolska-Kusnierz B, Pac M, Siewiera K, Galkowska E, Bernatowska E. B cell subsets in healthy children: Reference values for evaluation of B cell maturation process in peripheral blood. Cytometry Part B 2010; 78B: 372381.

Test developed and characteristics determined by ARUP Laboratories. See Compliance Statement A: aruplab.com/CS

H=High, L=Low, *=Abnormal, C=Critical

VERIFIED/REPORTED DATES

Procedure	Accession	Collected	Received	Verified/Reported
CD19+ B cells %	20-353-134185	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
CD19+ B cells	20-353-134185	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
CD20+ %	20-353-134185	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
CD20+	20-353-134185	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
Total Memory CD27+ %	20-353-134185	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
Total Memory CD27+	20-353-134185	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
Non switched CD27+IgD+IgM+ %	20-353-134185	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
Non switched CD27+IgD+IgM+	20-353-134185	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
Class-switched CD27+IgD-IgM- %	20-353-134185	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
Class-switched CD27+IgD-IgM-	20-353-134185	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
Transitional CD38+IgM+ %	20-353-134185	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
Transitional CD38+IgM+	20-353-134185	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
Plasmablasts CD38+IgM- %	20-353-134185	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
Plasmablasts CD38+IgM-	20-353-134185	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
Activated CD21low CD38- %	20-353-134185	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
Activated CD21low CD38-	20-353-134185	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00

END OF CHART

H=High, L=Low, *=Abnormal, C=Critical

Unless otherwise indicated, testing performed at:

ARUP LABORATORIES | 800-522-2787 | aruplab.com
500 Chipeta Way, Salt Lake City, UT 84108-1221
Tracy I. George, MD, Laboratory Director

Patient: Patient, Example
ARUP Accession: 20-353-134185
Patient Identifiers: 01234567890ABCD, 012345
Visit Number (FIN): 01234567890ABCD
Page 3 of 3 | Printed: 1/5/2021 9:04:23 AM
4848