

Client: Example Client ABC123
123 Test Drive
Salt Lake City, UT 84108
UNITED STATES

Physician: Doctor, Example

Patient: Patient, Example

DOB: 5/31/2007
Gender: Female
Patient Identifiers: 01234567890ABCD, 012345
Visit Number (FIN): 01234567890ABCD
Collection Date: 00/00/0000 00:00

Helicobacter pylori, Breath Test

ARUP test code 3006271

H. pylori, Breath Test

Positive * (Ref Interval: Negative)

INTERPRETIVE INFORMATION: Helicobacter pylori, Breath Test

A negative result does not rule out the possibility of H. pylori infection. If clinical signs are suggestive of H. pylori infection, retest with a new specimen or an alternate method.

Known causes of false-negative results include:
1. Ingestion of antimicrobials, proton pump inhibitors, and bismuth preparations during the preceding 2 weeks.

Known causes of false-positive results include:
1. Patients with achlorhydria.
2. Rinsing the testing solution in the mouth or not using the straw provided in the kit, which can allow contact with urease-positive bacteria.
3. The presence of other gastric spiral organisms such as Helicobacter heilmannii.

Note: The post-dose sample must be collected between 15 and 20 minutes post-dose to prevent a false-negative result.

VERIFIED/REPORTED DATES

Procedure	Accession	Collected	Received	Verified/Reported
H. pylori, Breath Test	23-098-102671	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00

END OF CHART

H=High, L=Low, *=Abnormal, C=Critical

Unless otherwise indicated, testing performed at: