

Client: Example Client ABC123
123 Test Drive
Salt Lake City, UT 84108
UNITED STATES

Physician: Doctor, Example

Patient: Patient, Example

DOB: 12/31/1996
Gender: Female
Patient Identifiers: 01234567890ABCD, 012345
Visit Number (FIN): 01234567890ABCD
Collection Date: 00/00/0000 00:00

QuantiFERON TB-Gold Plus, 1-Tube

ARUP test code 3017554

Quantiferon TB Gold Plus

Positive * (Ref Interval: Negative)

INTERPRETIVE INFORMATION:Quantiferon TB Gold Plus

Interferon gamma release is measured for specimens from each of the four collection tubes. A qualitative result (Negative, Positive, or Indeterminate) is based on interpretation of the four values: NIL, MITOGEN minus NIL (MITOGEN-NIL), TB1 minus NIL (TB1-NIL), and TB2 minus NIL (TB2-NIL). The NIL value represents nonspecific reactivity produced by the patient specimen. The MITOGEN-NIL value serves as the positive control for the patient specimen, demonstrating successful lymphocyte activity. The TB1-NIL tube specifically detects CD4+ lymphocyte reactivity, specifically stimulated by the TB1 antigens. The TB2-NIL tube detects both CD4+ and CD8+ lymphocyte reactivity, stimulated by TB2 antigens. An overall Negative result does not completely rule out TB infection.

A false-positive result in the absence of other clinical evidence of TB infection is not uncommon. Refer to: Updated Guidelines for Using Interferon Gamma Release Assays to Detect Mycobacterium tuberculosis Infection -- United States, 2010 (<http://www.cdc.gov/mmwr/preview/mmwrhtml/rr5905a1.htm>), for more information concerning test performance in low-prevalence populations and use in occupational screening.

Quantiferon Plus TB1 minus NIL **10.00 IU/mL H (Ref Interval: <=0.34)**

Quantiferon Plus TB2 minus NIL **10.00 IU/mL H (Ref Interval: <=0.34)**

QuantiFERON Mitogen minus NIL 10.00 IU/mL

QuantiFERON NIL 10.00 IU/mL

H=High, L=Low, *=Abnormal, C=Critical

Unless otherwise indicated, testing performed at:

VERIFIED/REPORTED DATES

Procedure	Accession	Collected	Received	Verified/Reported
Quantiferon TB Gold Plus	24-141-117326	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
Quantiferon Plus TB1 minus NIL	24-141-117326	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
Quantiferon Plus TB2 minus NIL	24-141-117326	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
QuantiFERON Mitogen minus NIL	24-141-117326	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
QuantiFERON NIL	24-141-117326	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00

END OF CHART

H=High, L=Low, *=Abnormal, C=Critical

Unless otherwise indicated, testing performed at:

ARUP LABORATORIES | 800-522-2787 | aruplab.com
500 Chipeta Way, Salt Lake City, UT 84108-1221
Jonathan R. Genzen, MD, PhD, Laboratory Director

Patient: Patient, Example
ARUP Accession: 24-141-117326
Patient Identifiers: 01234567890ABCD, 012345
Visit Number (FIN): 01234567890ABCD
Page 2 of 2 | Printed: 5/21/2024 5:58:08 PM
4848