

HOTLINE: Effective November 18, 2019

0010014 ABORH-PR **ABO-Rh Prenatal**

Performed: Mon-Fri Reported: 1-3 Days

 $\begin{array}{c} \textbf{Specimen Required: } \underline{Collect:} \ Lavender \ (K_2EDTA) \ or \ Pink \ (K_2EDTA). \\ \underline{Specimen \ Preparation:} \ Transport \ 7 \ mL \ whole \ blood. \ (Min: 3 \ mL) \end{array}$

Storage/Transport Temperature: Refrigerated.

Stability (collection to initiation of testing): Ambient: Unacceptable; Refrigerated: 1 week; Frozen: Unacceptable

 $\begin{tabular}{ll} \textbf{HOTLINE NOTE:} & \textbf{Remove information found in the Unacceptable Conditions field.} \end{tabular}$