

TEST CHANGE

Antibody Screen RBC with Reflex to Identification

0010020, ABSC-R

Specimen Requirements:

Patient Preparation:

Collect: Lavender (K2EDTA) or Pink (K2EDTA).

Specimen Preparation: Do not freeze. Transport 7 mL whole blood. (Min: 3 mL)

Transport Temperature: Refrigerated.

Unacceptable Conditions: Plasma Separator Tubes.

Remarks:

Stability: Ambient: Unacceptable; Refrigerated: 1 week; Frozen: Unacceptable

Methodology: Solid Phase Technology

Performed: Mon-Fri

Reported: 1-3 days

Note: If Antibody Screen is positive, Antibody Identification will be added. Additional charges apply.

CPT Codes: 86850; additional CPT codes may apply

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Reference Interval:

Test Number	Components	Reference Interval
	Antibody Screen Automated	Negative

HOTLINE NOTE: There is a reflexive pattern change associated with this test. One or more orderable or component has been added or removed to the reflexive pattern. Refer to the Hotline Test Mix for interface build information.