

HOTLINE: Effective November 18, 2019

0010020 Antibody Screen RBC with Reflex to Identification ABSC-R

Performed: Mon-Fri Reported: 1-3 Days

 $\textbf{Specimen Required:} \ \underline{\text{Collect:}} \ Lavender \ (K_2EDTA) \ or \ Pink \ (K_2EDTA).$

Specimen Preparation: **Do not freeze.** Transport 7 mL whole blood. (Min: 3 mL)

<u>Storage/Transport Temperature:</u> Refrigerated. <u>Unacceptable Conditions:</u> Plasma Separator Tubes.

Stability (collection to initiation of testing): Ambient: Unacceptable; Refrigerated: 1 week; Frozen: Unacceptable