
0020056**Viscosity, Serum****VIS-S**

Specimen Required: Collect: Serum separator or plain red tube.

Specimen Preparation: Transfer 1 mL serum to an ARUP Standard Transport Tube. (Min: 0.5 mL)

Storage/Transport Temperature: Refrigerated.

Unacceptable Conditions: Clotted specimens.

Stability (collection to initiation of testing): After separation from cells: Ambient: 24 hours; Refrigerated: 7 days; Frozen: 1 month

Reference Interval: Effective February 19, 2019

1.10-1.40 cP

Interpretive Data: See Compliance Statement B: www.aruplab.com/CS