

TEST CHANGE

Hepatitis Be Virus Antibody

0020095, HBEAB

Specimen Requirements:

Patient Preparation:

Collect: Serum separator tube (SST). Also acceptable: Lavender (EDTA) or green (lithium heparin).

Specimen Preparation: Separate serum or plasma from cells ASAP or within 2 hours of collection. Transfer 1 mL serum or plasma to an ARUP standard transport tube. (Min: 0.5 mL)

Transport Temperature: Refrigerated.

Unacceptable Conditions: Specimens that are heat-inactivated, grossly hemolyzed, grossly icteric, grossly lipemic specimens, or specimens containing particulate material.

Remarks:

Stability: After separation from cells: Ambient: ~~48~~24 hours; Refrigerated: ~~14~~7 days; Frozen: 30 days (avoid repeated freeze/thaw cycles)

Methodology: Qualitative Chemiluminescent Immunoassay (CLIA)

Note:

CPT Codes: 86707

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

This assay should not be used for blood donor screening, associated reentry protocols, or for screening human cells, tissues, and cellular- and tissue-based products (HCT/P).

Reference Interval:

Negative

Test Number	Components	Reference Interval
	Hepatitis Be Antibody	Negative