

Quarterly HOTLINE: Effective November 12, 2018

0020221 Hemoglobin, Urine HGBU

Specimen Required: Collect: Random urine.

Specimen Preparation: Centrifuge and separate urine from cells and other sediment. Transfer 4 mL aliquot of supernatant to an ARUP

Standard Transport Tube. (Min: 0.7 mL) Storage/Transport Temperature: Frozen.

Stability (collection to initiation of testing): Ambient: Unacceptable; Refrigerated: 8 hours; Frozen: 1 month