

TEST CHANGE

Hemosiderin, Urine

0020222, UHEMSID

Specimen Requirements:

Patient Preparation: First-morning collection is preferred.

Collect: Random urine.

Specimen Preparation: Mix specimen well. Transfer 4 mL to an ARUP ~~standard transport tube~~. **Standard Transport Tube**. (Min: 1 mL).

Transport Temperature: Frozen.

Unacceptable Conditions: Specimens in preservatives.

Remarks:

Stability: Ambient: 1 hour; Refrigerated: 1 week; Frozen: 1 week

Methodology: Semi-Quantitative Microscopy

Performed: Sun-Sat

Reported: ~~1-2 days~~
~~Within 24 hours~~

Note: Absent = Negative Present = 1+ to 4+

CPT Codes: 83070

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

As of Nov 16, 2020, reporting has switched from ~~semiquantitative~~ **semi-quantitative** to qualitative (absent versus present).

Reference Interval:

Effective November 16, 2020

Absent