

Quarterly HOTLINE: Effective February 19, 2019

0020414 Creatine Kinase Isoenzymes**CKISO**

Performed: Sun-Sat
Reported: 2-3 days

Specimen Required: Collect: Serum Separator Tube (SST) or Plain Red.

Specimen Preparation: Transfer 1 mL serum to an ARUP Standard Transport Tube. (Min: 0.5 mL)

Storage/Transport Temperature: Frozen.

Unacceptable Conditions: Specimens preserved in citrate, EDTA, fluoride, heparin, or iodoacetate. Grossly hemolyzed specimens.

Stability (collection to initiation of testing): After separation from cells: Ambient: Unacceptable; Refrigerated: 1 week; Frozen: 1 month (avoid repeated freeze/thaw cycles)

Note: This test will detect CK macroenzymes.