

HOTLINE: Effective February 22, 2022

0020850

Chloride, Urine

U CL

Specimen Required: Collect: 24-hour urine (without additives). Refrigerate during collection. Also acceptable: Random urine.

Specimen Preparation: Transfer 1 mL aliquot of urine from a well-mixed collection to an ARUP Standard Transport Tube. (Min: 0.5 mL) Record total volume and collection time interval on transport tube and test request form.

Storage/Transport Temperature: Refrigerated.

Stability (collection to initiation of testing): Ambient: 7 days; Refrigerated: 2 weeks; Frozen: 1 month