

TEST CHANGE

Factor VIII, Activity

0030095, F8

Specimen Requirements:

Patient Preparation:

Collect: Light blue (sodium citrate). ~~Special~~Refer to Specimen Collection and Handling Hemostasis/Thrombosis Specimens guide located at <https://www.aruplab.com/Specimen-Handling/SpecialSpecimenCollection/Hemostasis-Thrombosis.pdf> for hemostasis/thrombosis specimen handling guidelines.

Specimen Preparation: Transfer 2 mL platelet-poor plasma to an ARUP standard transport tube. (Min: 1 mL)

Transport Temperature: CRITICAL FROZEN. Separate specimens must be submitted when multiple tests are ordered.

Unacceptable Conditions: Serum. EDTA plasma, clotted or hemolyzed specimens.

Remarks:

Stability: Ambient: 4 hours; Refrigerated: Unacceptable; Frozen ~~at -20 Degrees C: 3 months; Frozen at -70 Degrees C: 6 months~~

Methodology: Electromagnetic Mechanical Clot Detection

Performed: Mon-Sat

Reported: 1-3 days

Note:

CPT Codes: 85240

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Reference Interval:

Age	Reference Interval
0-6 years	56-191%
7-9 years	76-199%
10-11 years	80-209%
12-13 years	72-198%
14-15 years	69-237%
16-17 years	63-221%
18 years and older	56-191%