

HOTLINE: Effective February 16, 2021

0050024

Albumin, Body Fluid

ALB-BF

Methodology: Quantitative **Immunoturbidimetry**/Quantitative Spectrophotometry

Specimen Required: Collect: Body fluid.

Specimen Preparation: Centrifuge and separate to remove cellular material. Transfer 1 mL body fluid to an ARUP Standard Transport Tube. (Min: 0.5 mL)

Storage/Transport Temperature: Refrigerated.

Remarks: Indicate source on the test request form.

Unacceptable Conditions: Contaminated or grossly hemolyzed specimens. Needle sent with specimen.

Stability (collection to initiation of testing): After separation from cellular material: Ambient: **Unacceptable**; Refrigerated: **1 month**;
Frozen: 1 month (avoid repeated freeze/thaw cycles)

HOTLINE NOTE: There is a clinically significant charting name change associated with this test.

Change the charting name for component 0050024, Albumin Level Body Fluid from Albumin Level Body Fluid to **Albumin, Body Fluid**.