
0050049 **Immunofixation Electrophoresis, Immunoglobulin D and Immunoglobulin E, Serum** **IFE D/E**

Performed: Sun-Sat
Reported: 1-5 days

Specimen Required: Collect: Serum separator tube.

Specimen Preparation: Separate serum from cells ASAP or within 2 hours of collection. Transfer 1 mL serum to an ARUP Standard Transport Tube. (Min: 0.5 mL)

Storage/Transport Temperature: Refrigerated.

Unacceptable Conditions: Plasma.

Stability (collection to initiation of testing): After separation from cells: Ambient: **Unacceptable**; Refrigerated: 1 week; Frozen: 1 month